

# Call cord stations with reset button

# Traditional and dependable solution for bedside emergency calls

Lifeline 900 MHz wireless nurse call cords are designed to be used with EchoStream® and Inovonics backbones for bedside emergency calls and patient requests. Used in a community setting, help requests from connected call cords are processed by the CarePoint resident safety system, which receives alerts and notifies supervising staff — helping you build a better community experience.

#### Features and benefits

Local reset button allows staff to acknowledge alarm calls at source

Removable 10-foot call cord features momentary push button

**Supervised wireless transmitter** provides notifications for service needs or battery replacement

LED Indicator on enclosure that illuminates for 10 seconds when activated and reset

Easy to clean plastic pull cord and membrane switch

Long-life 3.6 volt lithium battery delivers up to five years of battery life

**UL 2560 certified** 

### Specifications

Color	White with grey switch
Dimensions	5" H x 3.125" W x 1.7" D
Weight	13.1 oz
Battery	Panasonic CR 123A 3V lithium
Battery life	Up to 5 years
Check-in time frequency	60 minutes
Frequencies	EchoStream®: 902-928 MHz, frequency hopping spread spectrum
	Lifeline: 900 MHz
Regulatory compliance	FCC Part 15, Industry Canada, RoHS, UL 2560
Operating environment temperature	-20° to 140°F (-20° to 60°F)
Humidity	Up to 90% (non-condensing)

<sup>\*</sup>Performance is a function of many variables, including repeater layout, device type and quantity, as well as the structure and environment in which the network operates. Continual enhancements to products may cause specifications to change without notice. Patent pending.

## Ordering information

Description	Part number
900MHz non-locking cord station with reset	300005984471
EN1341CC Call cord station, with reset	300005399101
EN1341CC Call cord station, no reset	1902060
EN1341CC Call cord station, no reset, with non-locking cord	1902061
EN1341CC Call cord station with reset, with non-locking cord	1902063

To order parts and accessories, please contact Lifeline Senior Living at **1-800-816-4885**, option 2.



1-800-387-8120

Toronto, Ontario, M3C 3E9 Canada