

Personal help buttons

Encourage constant wear for continuous protection

Lifeline 319MHz and 900MHz personal help buttons are lightweight and waterproof personal help devices. They can be worn by residents as a pendant with an adjustable neck cord or as a wristband with a soft, elastic strap. As an essential part of CarePoint resident safety systems, these help buttons provide maximum comfort and flexibility to encourage consistent wear by residents – helping you build a better community experience.

Features and benefits

Comfortable, waterproof design encourages constant wear and continuous protection

Bright LED indicates that a help call has been sent

Special button design reduces accidental activations

Supervised device provides notification for servicing or low-battery condition

Elastic wristband or adjustable lanyard with breakaway feature helps to improve comfort and safety

Specifications

Frequency	319.5MHz / 900MHz	
Weight	0.9 oz (25.7 g)	
Materials	Button: Silicone	
	Lanyard: Nylon	
	Wristband: Thermoplastic elastomer	
Indicator	Green LED	
Signal range	600 ft (182.8 m) open field range	
Battery life	Five years*	
Product limited warranty	Two years	
Operating environment	Temperature: 23 to 122°F (-5 to 50°C)	
	Waterproof to IP67 standard	
Regulatory compliance	FCC, RoHS, IP67, conforms to IC standards	

Note: The range and performance of any wireless product depends on the structure and environment in which it operates. Continual enhancements to Lifeline products may cause specifications to change without notice.

*Battery life is calculated at an average of one call placed per day.

Ordering information

Description	Part number
Pendant with neck lanyard and wristband (319PHB), black	1902023
Pendant with neck lanyard and wristband (319PHB), white	1902025
Pendant with neck lanyard and wristband (900PHB), black	1902034
Pendant with neck lanyard and wristband (900PHB), white	1902035
Replacement wrist band, black	1902043
Replacement wrist band, white	1902044
Replacement lanyard, black	1902045



To order parts and accessories, please contact Lifeline Senior Living at **1-800-816-4885**, option 2.

Lifeline

200 Donald Lynch Blvd, STE 300 Marlborough, MA 01752 1-800-816-4885

Lifeline Canada

95 Barber Greene Road, Suite 105 Toronto, Ontario, M3C 3E9 Canada 1-800-387-8120 © 2022 Lifeline. All rights reserved. Trademarks are the property of their respective owners. www.lifeline.com/

seniorliving MKT SL007 REV3