

Personal help buttons with fall detection

Smart pendant-style buttons that automatically send an alarm if a fall is detected

Lifeline 319MHz and 900MHz fall detection pendants are lightweight and water-resistant, encouraging constant wear and continuous protection. This vital part of the CarePoint resident safety system alerts staff to potential falls – even if a resident is unable to press their help button. And it's one more way you can count on Lifeline to help you build a better community experience.

Features and benefits

Lightweight, water-resistant design encourages constant wear and continuous protection

Bright LED indicates that a help call has been sent

Special button design reduces accidental activations

Adjustable lanyard with breakaway feature enhances user safety

Supervised device provides notification for low-battery condition

Specifications

| Frequency | 319.5MHz / 900MHz | |
|--------------------------|--|--|
| Weight | 1.01 oz (28.6 g) | |
| Materials | Button: Silicone | |
| | Lanyard: Nylon | |
| | Pendant Enclosure: ABS plastic | |
| Indicator | Green LED | |
| Low battery indicator | Red LED | |
| Signal range | 600 ft (182.8 m) open field range | |
| Battery type | 3V Lithium (CR2477), replaceable | |
| Battery life | 1.5 years* | |
| Product limited warranty | Two years | |
| Operating environment | Temperature: 23 to 122°F (-5 to 50°C) | |
| | Liquid and Solid Ingress Protection Rating: IP45 | |
| Regulatory compliance | FCC, RoHS, IP45, conforms to IC standards | |

Note: The range and performance of any wireless product depends on the structure and environment in which it operates. Continual enhancements to Lifeline products may cause specifications to change without notice.

Lifeline Senior Living technology does not detect 100% of falls. If able, users should always push their button when they need help.

*Battery life is calculated at an average of one call placed per day.

Ordering information

| Description | Part number |
|------------------------------|-------------|
| Fall detection pendant (319) | 1902024 |
| Fall detection pendant (900) | 1902036 |
| 1902036 Replacement lanyard | 1902045 |

To order parts and accessories, please contact Lifeline Senior Living at **1-800-816-4885**, option 2.

Lifeline

200 Donald Lynch Blvd, STE 300 Marlborough, MA 01752 1-800-816-4885 Lifeline Canada 95 Barber Greene Ro

95 Barber Greene Road, Suite 105 Toronto, Ontario, M3C 3E9 Canada 1-800-387-8120 © 2022 Lifeline. All rights reserved. Trademarks are the property of their respective owners. www.lifeline.com/seniorliving