

PHILIPS

Carepoint

Call Cord

with Reset Button



Wireless Nurse Call button with local reset for closed loop response. Call cords are available with 7 or 10 foot cords and are available either for flush or surface mount.

Features & Benefits

Traditional and dependable solution for bedside Emergency Call

Available in Flush or Surface mount models for ease of installation in new and retrofit projects

Self contained unit with Supervised Wireless transmitter

Long-Life 3.6 Volt Lithium battery, changeable on-site

Unit Reset Button for acknowledging alarm calls at source

Removable 7 ft. or 10 ft. call cord with momentary push button

Supervised device which provides notification when device needs service or battery replacement

Specifications

Color	White face plate with black button reset switch
Dimensions	Flush Mount: 4.5" L x 2.75" W x 1.125" D, with cord installed Surface Mount: 4.75" L x 3.125" W x 2.125" D, with cord installed
Call Cord	7 ft. or 10 ft, gray color with momentary push button, 1/4" phone plug & security clip
Battery	3V DC Lithium
Average Battery Life	5 years

Mounting Note: Flush Mount devices require a non-metallic single 18.0 cu. in. gang box

